



**PCI**  
EUROPE

## **INSTRUCTIONS FOR WATER ANALYSIS TEST**

### WATER ANALYSIS TEST

Please read instructions completely before you begin.

Unscrew cap to container provided.

1. Provide approximately 70ML of water specimen in container (if water is standing), sample from deepest possible safe area beneath the surface.
2. Screw cover back on container.
3. Fill out your name and location of where water sample was taken from (under donor's initials).
4. Store in a cool place prior to mailing back to PCI Europe for testing.
5. Results will be made available within 7-10 days of receipt of the water sample.

**Please send test kit back in original packaging**

If you have any questions please contact us on [enquiries@parasitetesting.co.uk](mailto:enquiries@parasitetesting.co.uk).

Broadlands Farm  
101 Durkar Low Lane  
Durkar  
Wakefield  
WF4 3BJ  
01924 242851